


PORKKA



Installation and operation manual

Vaccine Storage Rooms



Vaccine Storage Rooms Installation and Operation Manual

TABLE OF CONTENTS

1.	GENERAL	3
1.1.	Placement (ventilation)	3
2.	SAFETY	4
2.1.	Occupational safety	4
3.	RECEIVING AND STORING THE DELIVERY	6
3.1.	Receiving the delivery	6
3.2.	Storage	6
3.3.	Disposal of packaging material	7
4.	VACCINE STORAGE ROOM MODELS	8
4.1.	Delivery content and main parts	8
5.	MOUNTING	10
5.1.	Room assembly	10
5.2.	Mounting cold room without ceiling	11
5.3.	Mounting refrigeration cooling units	11
5.4.	Mounting ceiling	12
5.5.	Mounting control box unit	12
5.6.	Mounting Alarm Unit	13
5.7.	Mounting temperature sensor	13
5.8.	Mounting Dickson KT 7 temperature chart recorder (optional)	14
5.9.	Routing of cables	14
5.10.	Mounting lamps	15
5.11.	Installing/connecting electricity system	16
5.12.	1-phase wiring diagram	17
5.13.	3-phase wiring diagram	18
5.14.	3-phase wiring diagram automatic duty sharing with a timer relay	19
5.15.	Alarm Unit Lite 2018 Battery Backup	20
5.16.	Lay-out of control box electricity	21
5.17.	Installing/connecting mains power supply (tästä valokuva)	22
5.18.	Installing shelves	22
5.19.	Testing and commissioning	23
6.	OPERATION	23
6.1.	Placing items in the cold/freezer rooms	23
6.2.	Operation	23
7.	FEATURES	24
7.1.	Electrical connections	24
7.2.	Operational faults	24
7.3.	Combined electronic control board and temperature display	25
7.4.	Keys	25
7.5.	Function of key lamp lights	26
7.6.	Function of the display lights	26
7.7.	Stand-by	26
7.8.	Adjusting the temperature	26
7.9.	Defrost	26
7.10.	Manual defrost	26
7.11.	Checking minimum temperature	27
7.12.	Checking maximum temperature	27
7.13.	Resetting maximum and minimum temperature recorded	27
7.14.	Alarm signals	27
8.	MAINTENANCE	28
8.1.	Cleaning	28
8.2.	Guarantee	29
8.3.	Disposing of the unit	29
9.	LAY-OUT OF WIC 10	30
10.	LAY-OUT OF WIC 20	31
11.	LAY-OUT OF WIC 30	32
12.	LAY-OUT OF WIC 40	33
13.	LAY-OUT OF WIC/WIF 40	34

Vaccine Storage Rooms Installation and Operation Manual

1. GENERAL

These instructions are intended for installation and operation of vaccine storage cold and freezer rooms.

It is very important to read these installation and operating instructions carefully before you start using the equipment for the first time. Also, please keep these instructions in a safe place for future reference or use by another operator.

Modular cold and freezer vaccine storage rooms are manufactured for the storage of chilled and frozen items. They are not designed for freezing or chilling down of hot foodstuffs, or other items.

In this instruction manual, measures have been described relating to the installation and daily use of the product. Yearly maintenance and repairs have to be done by an authorized service company. All the covers of the product have to be properly fastened in place before switching on the product.

By following these instructions you can improve the product's performance and reduce unnecessary repair costs. Please note: it is important to have your equipment regularly maintained by a professional engineer.

You will find the terms of guarantee in this manual (see the Table of Contents).

1.1. Placement (ventilation)

Before erecting the room, ensure that there is sufficient ventilation in the area where the room is to be located. During operation, the room will dissipate a heat load of 1 - 2,1 kW per hour in the area in question. This heat must be removed from the area by means of ventilation, such as louvers or extraction systems. If impractical, consult your supplier for the possibilities of using our heat disposal system. The refrigeration systems are designed to operate satisfactorily in ambient max. +43 °C. Ensure that the floor area to be used is level. Requirement of the evenness is max ± 3 mm/m. Correct any unevenness at this time. It is recommended to leave a gap of 50 mm minimum between a freezer room and the surfaces of the surrounding building for air circulation purposes.

The floor under the freezer room should be treated against moisture and ice. If needed, the floor elements can be delivered with floor heating elements (optional accessory). Before unpacking, ensure that there is sufficient space where the room is to be assembled, including sufficient height.

Vaccine Storage Rooms Installation and Operation Manual

2. SAFETY

DANGER!

Follow general rules and regulations concerning occupational safety.



2.1. Occupational safety

Always wear a safety helmet, safety shoes, protective gloves, and protective clothing in the installation work.



Be careful when working near lifting devices.



Vaccine Storage Rooms Installation and Operation Manual

Before use, ensure proper operation and condition of scaffolds, ladders, and lifting devices and equipment.




Caution!

Risk of injury! Be extra careful during lifting work.




Electrical appliances are to be installed and repaired only by experts qualified in maintaining electrical devices.



Vaccine Storage Rooms Installation and Operation Manual

3. RECEIVING AND STORING THE DELIVERY


3.1. Receiving the delivery



Attention!

The manufacturer's guarantee does not cover damages caused by transportation!

3.2. Storage





Vaccine Storage Rooms Installation and Operation Manual

3.3. Disposal of packaging material


Note!

The packaging material must be recycled as efficiently as possible given the resources of the installation site.



Vaccine Storage Rooms Installation and Operation Manual

4. VACCINE STORAGE ROOM MODELS



4.1. Delivery content and main parts

The delivery for a walk-in cold and freezer room includes:

Wall, ceiling, and floor panels

Panel mounting accessories

Pressure-balancing valve

- Provided as standard for freezer rooms and by special order for cold rooms as well

Parts list of elements, door installation and installation materials for panels (See Panel Unijoint manual)

Doors

Floor plan

Refrigeration units

Electrical accessories

Temperature logging accessories (Temperature recorder) (optional)

Special drawings

Mounting instructions

Lights

Spare parts of Euromon

Control box

Special service equipment

Refrigerant, refrigerant units

Shelvings

Valve set (Only in WIC WIF 40)

Vaccine Storage Rooms Installation and Operation Manual

Example of a delivery parts list.

PARTS LIST					
Qty.	Code	Size	Description	Note	Crate
VALVE SET (FACTORY INSTALLED IN PANEL W1)					
1	PC	AP75	CONNECTION BOX (FOR PRESSURE REL.VALVE) KYTKENTÄRASIA		} 1
1	PC	BQ18	PTV MINI ELEBAR LT 230V 50HZ 7w MINIELEBAR PRESSURE RELIEF VALVE, 230 V		
6	PC	AM69	LEVYRUUVI LIERIÖK. 4,2x16 A4 PLATE SCREW 4,2x16 A4	DIN 7981 RISTIÜ	
INSTALLATION MATERIAL FOR PANELS					
70	PC	AN76	PEITET.9,5-11 TUCKER P-208/0 VALK PLASTIC PLUG,WHITE		} 5
10	PC	AS86	BUTY-MASSA VALKOINEN PLASTINEN BUTYL SEALANT WHITE 0.31 L		
1	PC	AS87	SEALANT PRESS 1/3 L SEALANT PRESS 1/3 L		
2	PC	AS99	SODAL SILIRUB N 0.3L VÄRITÖN Silicon neutral		
10	m	AY97	TIIVISTE PE-PROF.20X68 A29702/c-2 GASKET, SOFT 20x68 mm 7 ST L=2400 mm		
15	m	BA21	TIIVISTE KOVA PU-50, 20X68 MM GASKET,HARD 20X68 mm 8 ST L=2000 mm		
2	PC	BC49	PIKALUKON AVAIN 100 MM A33979/a CAM-LOCK KEY 100 mm 30445		
2	PC	CA35	TREMSIL 500 0.3L valkoinen TREMSIL 500 0.3L WHITE SILICON E SEALANT		
2	PC	CD69	TREMSIL 500 0.3L HARMAA TREMSIL 500 0.3L GREY SILICONE SEALANT		
1	PC	YY90	PIIRUSTUS JA OSALUETTELO WALK-IN MOUNTING AND OPERATION MANUAL		
1	PC	UN01	LIGHTS		
15	M	AK37	WIRE MMJ 2X1,5 N 04064127 SÄHKÖJOHTO		} 5
2	PC	AE85	LIGHT EQUIPMENT AVR 10 E27/60 W VALAISIN	4211418	
2	PC	AJ74	BULB 60 W H 230V E27 VAKIOLAMPPU	4731413	
1	PC	BT78	LIGHT SWITCH 2-N MV ORANGE VIPUPAINOKYTKIN	2004282	
1	PC	BT81	PASS-THROUGH PIPE 16x95 A13684-6 LÄPIVIENIPIUTKI		
2	PCS	AJ37	PIPE END COVER JAPP-16, JPP-16 PUTKENPÄÄTE	1131013	
50	PCS	AX83	FASTENER FOR WIRE TK10-14 JOHDON KIINNITIN	1310114	
50	PCS	PR60	SCREW 4.2X13 SZN WRONIC-RUUVI		
8	PCS	AM89	PLATE SCREW 4.2X22 SZN LEVYRUUVI ITSEPÖR.	DIN 7504-N, DIN 7981	

38470


Drawn	TTa	Order	550933	Proj.m:	TOM	Room Number:	550933-01
Drawn date	2015-12-07	Project					
Delivery from factory	2016-01-14	Room	45143637, Customer Ref.:				
PORKKA PORKKA FINLAND OY P.O. Box 127 FIN-33101 Tampere Finland Puh. 020 5555 12 Tel. +358 20 5555 12 Fax +358 20 5555 360 www.huurre.com							
		Drawing number 550933-01 -05					2/5

Vaccine Storage Rooms Installation and Operation Manual


5. MOUNTING

5.1. Room assembly


Assemble the room according to the following assembly instructions.




PANELS UNIJOINT




HMSS SHELVES



S01-DOOR



REFRIGERATION UNIT



DICKSON KT (optional)
See temperature recorder at
www.dicksondata.com



GLAMOX I25-800
LED 2400 HF 840 OP
Light

Vaccine Storage Rooms Installation and Operation Manual



5.2. Mounting cold room without ceiling

Attention!

Make sure the floor is flat/even!


Note!

See the list of 'Required equipment' and 'Most common mounting accessories' in document WI_Rooms_Manual_Panels_Unijoint_EN.pdf




5.3. Mounting refrigeration cooling units

Please see the refrigeration Installation Instructions.



Vaccine Storage Rooms Installation and Operation Manual

5.6. Mounting Alarm Unit



Part No.	Description	Qty	Unit	Notes	Drawn
PARTS LIST					
FACTORY INSTALLED (IN PANEL KIT)					
1	ALARM UNIT	1	UNIT	1	
2	ALARM BATTERY	1	BATTERY	1	
3	ALARM BATTERY COVER	1	COVER	1	
4	ALARM BATTERY COVER SCREW	2	SCREW	2	
5	ALARM BATTERY COVER GASKET	1	GASKET	1	
6	ALARM BATTERY COVER O-RING	1	O-RING	1	
7	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
8	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
9	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
10	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
11	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
12	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
13	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
14	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
15	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
16	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
17	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
18	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
19	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
20	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
21	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
22	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
23	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
24	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
25	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
26	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
27	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
28	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
29	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
30	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
31	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
32	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
33	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
34	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
35	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
36	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
37	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
38	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
39	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
40	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
41	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
42	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
43	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
44	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
45	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
46	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
47	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
48	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
49	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
50	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
51	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
52	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
53	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
54	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
55	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
56	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
57	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
58	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
59	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
60	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
61	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
62	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
63	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
64	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
65	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
66	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
67	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
68	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
69	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
70	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
71	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
72	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
73	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
74	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
75	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
76	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
77	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
78	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
79	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
80	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
81	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
82	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
83	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
84	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
85	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
86	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
87	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
88	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
89	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
90	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
91	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
92	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
93	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
94	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
95	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
96	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
97	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
98	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
99	ALARM BATTERY COVER GROMMET	1	GROMMET	1	
100	ALARM BATTERY COVER GROMMET	1	GROMMET	1	

5.7. Mounting temperature sensor

Temperature sensor is included in alarm unit.
Install the sensor inside the cold room and the location must be in the center area vertically.

Note!

Sensor should not be behind the shelving or near the evaporator, or door opening.



Vaccine Storage Rooms Installation and Operation Manual

5.10. Mounting lamps

Only a qualified electrician is allowed to open the luminaire dome for installation and adjustment purposes. In case of breakage, only the manufacturer is allowed to repair the product.

It is recommended to install the cables for the lamp on the outside surface of the room. The cables can also be installed inside the room, if needed (see the following page).


For further information on mounting and adjusting the lamp, see light fitting mounting instructions.



PARTS LIST				
Type	Description	Qty	Unit	Order
1	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
2	POWER CABLE FOR EXTERNAL CONNECTION	1		1
3	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
4	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
5	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
6	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
7	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
8	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
9	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
10	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
11	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
12	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
13	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
14	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
15	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
16	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
17	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
18	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
19	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1
20	CONNECTION KIT FOR EXTERNAL CONNECTION	1		1


Vaccine Storage Rooms Installation and Operation Manual

5.11. Installing/connecting electricity system



Vaccine Storage Rooms Installation and Operation Manual


5.12.1-phase wiring diagram



Vaccine Storage Rooms Installation and Operation Manual


5.13. 3-phase wiring diagram

30926_V2



Vaccine Storage Rooms Installation and Operation Manual


5.14.3-phase wiring diagram automatic duty sharing with a timer relay



Vaccine Storage Rooms Installation and Operation Manual

5.15. Alarm Unit Lite 2018 Battery Backup

30929_V10_2



NORMAL SITUATION:

POWER SUPPLY OK, ELECTRICITY SWITCHED ON,
DOOR CLOSED, TEMPERATURE INSIDE LIMIT VALUE

ALARM / MALFUNCTION:


TIME RELAY (KT1) WILL START COUNTING:
IF ELECTRICITY IS SWITCHED OFF, OR
IF POWER SUPPLY IS MISSING,
OR IF THERMOSTAT IS NOT OK
TIME RELAY (KT2) WILL START COUNTING:
IF DOOR IS OPEN

WHEN TIME RELAY (KT1/KT2) REACHES THE ADJUSTED TIME,
TIME RELAYS (KT1/KT2) CONTACT CHANGES POSITION AND ALARM IS ON

ALARM CAN BE RESET WITH THE RESET SWITCH. IF THE CAUSE OF THE ALARM HAS NOT BEEN
REMOVED,
TIME RELAY WILL GIVE THE ALARM AGAIN AFTER THE ADJUSTED TIME

IF ELECTRICITY IS SWITCHED OFF OR POWER SUPPLY FAILS,
BATTERY T1.2 12Ah PROVIDES POWER FOR 1 DAY

30929_V10_3





+ 43 °C
IEC-60335-1

Vaccine Storage Rooms Installation and Operation Manual


5.16. Lay-out of control box electricity

30928_V7




Vaccine Storage Rooms Installation and Operation Manual

5.17. Installing/connecting mains power supply



5.18. Installing shelves

Please see Manual_WI-rooms_Shelving_HMSS.pdf



Vaccine Storage Rooms Installation and Operation Manual

5.19. Testing and comissioning

1. Ensure that the power cable is connected to the mains
2. Turn the main switch to ON position/Turn the main switch on
3. Set Control switch to 1 on the control unit
4. Ensure that all motorised fans rotate freely,
5. Set the temperature at +4 °C for the positive unit and -20 °C for the negative unit,
6. Press the ON/OFF button on the refrigeration unit regulator
7. Refrigeration unit No. 1 starts up when it must be colder,
8. Emergency stop: Ensure that the main power switch button is working properly,
9. Check that the lighting inside the cold room functions properly.
10. Set the Manual Change over switch to position No. 2 to test that refrigeration unit No. 2 functions correctly (in the event that refrigeration unit No. 1 is malfunctioning or requires maintenance),
11. Press the ON/OFF button on the refrigeration unit regulator,
12. Repeat the measures you made to refrigeration unit No. 1 to refrigeration unit No. 2.
13. Set the temperature to +4 °C for the positive unit and -20 °C for the negative unit.

6. OPERATION

6.1. Placing items in the cold/freezer rooms

When placing items, especially in the freezer room, pay special attention to sufficient free air space around each package to ensure proper air circulation. Do not place or stack packages or boxes in front of the air movement fan. Items should be stored in closed packages.

6.2. Operation

All refrigeration units are pre-tested in the factory and the operational temperatures can be adjusted as follows:

- | | | |
|----------------------|--------------|-----------------|
| a) cold rooms (C) | +2...+8 °C | (+36... +46 °F) |
| c) freezer rooms (F) | -28...-15 °C | (-18... +5 °F) |

Before using the room, check that the operational temperature is reached. If the operational temperature is not reached, see the "Operational Faults" section of these instructions before you call in a maintenance engineer. A strong smell caused by cleaning agent or silicon mastic may appear in the room, in which case the room should be washed with a mild detergent suitable for the purpose, then dried and ventilated for at least 24 hours.

Do not remove any panels or guards or the electrical control panels before disconnecting the unit from the mains supply by removing the plug or the fuse from the socket!

Vaccine Storage Rooms Installation and Operation Manual

7. FEATURES

7.1. Electrical connections

Before connecting the machinery to the mains electricity, ensure that the voltage is the same as that shown on the refrigeration unit serial plate. Ensure that a separate supply is provided for each room being installed. Refrigeration unit is connected to a control box unit.

Please see manual of refrigeration unit's technical data for power supply and fuse.

The refrigeration unit is supplied with a flexible cable with a molded electrical plug top with fuse. Three phase models require an isolating switch fuse.

In conformity with EN 61000-3-11, the maximum value of supply impedance Z_{max} is 0,27ohm at Flicker test.

7.2. Operational faults

If the operational temperature is not reached or the alarm is activated, check the following:

- The door has not been left open for a long period.
- Electrical supply has not been cut off.
- Defrosting is not switched on.
- The equipment is not overloaded with hot items.
- An attempt has not been made to get the equipment colder than the manufacturer's stated operational temperature.
- Excess ice has not formed on the evaporator. If this has happened, carry out a complete defrosting.
- The "CSD"-warning is not displayed.
- The ambient temperature is not too high nor too low.
- The unit is not in STAND-BY -mode ("OFF" on display).

If the operational faults continue after you have checked the above-mentioned points, transfer the stored items to an appropriate back-up storage to prevent them from spoiling and call your refrigeration supplier or service agent.

Vaccine Storage Rooms Installation and Operation Manual

7.3. Combined electronic control board and temperature display

The equipment is fitted with a combined digital temperature display with many built-in features, including data connections for monitoring of equipment temperatures and alarms. The potential-free relay is either normally open or closed for a connection to a building maintenance alarm system or a supervising system.

During normal operation, the display indicates the internal air temperature of the room. When 'DEF is displayed, the room's automatic defrost period is ongoing.

7.4. Keys



The key shows the highest thermostat value registered. In the programming mode, you navigate inside the program or use the key to raise the selected value.



The key shows the lowest thermostat value registered. In the programming mode, you navigate inside the program or use the key to lower the selected value.



Manual defrosting starts when you press and hold the key for 3 seconds.



The key for checking and changing the settings. In programming mode, the key is used to select the value and to check the programming function. The registered value of the highest and lowest temperature can be deleted by pressing the key for 3 seconds when the value is displayed.



Light switch (Note! Used only in cabinets with lights)







Power switch










The keylock prevents unintended use of the keys. The lock is switched on by holding the keys down at the same time for 3 seconds. The display momentarily shows the letters "PoF". With the keylock on, you can still check the lowest and highest temperature registered. The light switch function can still be used when the keypad is locked (only on cabinets fitted with lights). To unlock the keyboard, press the keys at the same time for more than 3 seconds. The "Pon" message will be displayed and the keyboard is unlocked.

Vaccine Storage Rooms Installation and Operation Manual



7.5. Function of key lamp lights

KEY	MODE	FUNCTION
	Green light on	Manual defrosting is on
	Green light flashing	Set temperature is displayed and can be changed
	Green light on	Light inside the cabinet is on
	Red light on	The equipment is on standby mode, ON/OFF





7.6. Function of the display lights

LIGHT	MODE	FUNCTION
	On	Compressor is running
	Flashing	Time delay to protect compressor after start-up
	On	The fan of the evaporator is on
	On	Defrosting is on
	Flashing	Time delay after defrost
	On	ALARM
	On	The fan of the condenser is on

7.7. Stand-by

When pressing the  (ON/OFF) button the controller will display 'OFF' for 5 seconds. After this, the light at the top right hand of the switch will be illuminated. Note: when switched 'OFF', all relays are shut down and the unit will not record either the controller data or alarms. Note: The light switch will still operate as normal. To switch the power on (for example, after deep clean), simply press the  ON/OFF button again.


7.8. Adjusting the temperature

Push and immediately release the  key. The display will show the set temperature. You can then change the set temperature by pushing either the  or  key to the temperature required. Note: this must be done within 10 seconds. To memorise the new set temperature, press the  key again. After this the newly set temperature and the set key light will start to flash, or wait for 10 seconds when the display reverts to the true temperature within the equipment.

7.9. Defrost



The equipment is designed to initiate automatic defrosts. During this period, 'DEF' is displayed on the controller. After the defrost period, the unit will revert to display the actual temperature. Chiller units defrost every 12 hours, and freezer versions defrost every 6 hours.

7.10. Manual defrost



Press the  key for more than 2 seconds to start a manual defrost. Manual defrost resets the automatic defrost counter, and the unit continues to operate automatically after the defrost has completed. Initiating a manual defrost may be necessary for instance when the door has been left open by accident, causing extra formation of ice.

Vaccine Storage Rooms Installation and Operation Manual

7.11. Checking minimum temperature


Press and release the  key. "Lo" will appear on the display, followed by the minimum temperature logged. Press the  key again or wait for 5 seconds and the display will revert to the actual temperature within the equipment.

7.12. Checking maximum temperature

Press and release the  key. "Hi" will appear on the display, followed by the highest temperature logged. Press the  key again or wait for 5 seconds and the display will revert to the actual temperature within the equipment.

NOTE! After installation, reset the temperature stored.

7.13. Resetting maximum and minimum temperature recorded

To reset the stored temperature, when maximum or minimum temperature is displayed: Press the  key until the "rST" label starts blinking.

7.14. Alarm signals

CODE	REASON
"HA"	HIGH TEMPERATURE ALARM The equipment has detected high temperature (warm items have been placed in the equipment, door left open, etc.) The alarm will automatically stop when the temperature reverts to the set temperature, or defrost is initiated.
"LA"	LOW TEMPERATURE ALARM Check that the items stored within the equipment are at the required temperature. If necessary, move the stored items to safety to prevent any damage to them. Initiate manual defrost (see chapter Manual defrost, page 21). If the alarm continues, silence the buzzer relay and contact your authorised service dealer immediately.
"cSd"	CONDENSER TEMPERATURE ALARM The alarm is caused by high condensing temperature (e.g., abnormally high ambient temperature or condenser blocked). Put the controller on standby and unplug or disconnect the equipment from the electrical supply. Clean the condenser filter by washing or replacing if required. Also clean the finned condenser with a soft brush or vacuum as necessary.
"EE"	DATA OR MEMORY FAILURE * The controllers are provided with an internal checking system for data integrity. The alarm 'EE' flashes when a failure in memory data occurs. In this instance the alarm output is activated. To reset 'EE' alarm status and restart normal function, press any button. The display will show 'rSt' for approximately 3 seconds.

* The control unit will operate in backup mode. Reset the alarm and contact the authorised service dealer immediately.

NOTE! If necessary, remove stored items to a back-up storage as required. Porkka and/or its distributors and/or its service contractors will not, under any circumstances, be liable for deterioration or loss of stored items, however it may occur. Stored items should be insured by the equipment operator.

Vaccine Storage Rooms Installation and Operation Manual

8. MAINTENANCE

8.1. Cleaning

Do not use detergents or disinfectants containing chlorine, solvents, scrubbing products, a knife, or other sharp tools.

Refer to the product description of the disinfectant to see which materials it is suitable for. Do not let it splash on sensitive parts, such as thermostat regulator or door hinges. Dry the equipment after disinfecting and leave it to ventilate.

Room with standard floor:

Vaccine storage rooms should be completely defrosted two to three times a year. Clean all surfaces in the room (even the shelves) carefully with a wet cloth and a mild water soluble detergent. Do not use corrosive detergents including chlorine or acetic acid.

Remove detergent remnants by cleaning the surfaces once again with a wet cloth and clean water. Wipe the surfaces dry with a dust-free cloth. Do not use running water or a power washer to clean the room.

In connection with this comprehensive cleaning, it is useful to organize the periodic overhaul of the refrigeration equipment. Regular yearly maintenance guarantees long trouble-free operation of the equipment. Maintenance carried out periodically will save energy and reduce running costs. Inspecting correct operation of the equipment is the user's responsibility.

Room without standard floor:

Rooms without standard floor should be defrosted and serviced according to the instructions above. Running water can be used to clean the floor. After the floor has been cleaned, wipe all surfaces dry with a dust-free cloth. Do not use a pressure washer or corrosive detergents.

Put back all parts removed during cleaning and switch on the equipment.

Ensure that the temperature has returned to normal storage level before using the equipment again.

Vaccine Storage Rooms Installation and Operation Manual

8.2. Guarantee

Check the guarantee period from your supplier.

The guarantee does not cover faults caused by:

- Transportation
- Overloading or user's negligence
- Negligence due to not reading manuals, negligence of proper care and maintenance
- Changes in current (max $\pm 10\%$ allowed) caused for example by lightning etc.
- Modifications or repairs performed by an unauthorized service agent
- Use of parts not supplied and approved by the manufacturer

The guarantee does not cover:


- Incidental scratches/marks or other minor faults which appear during unpacking or installation and which do not affect the operation or performance of the equipment.

The manufacturer or their selling agent is not at any time, or under any circumstances, liable for the loss of stored items. The owner/user of the equipment must ensure that the contents are insured at all times.

All goods are supplied under our terms and conditions of sale, a copy of which may be obtained upon request.


8.3. Disposing of the unit

Dispose of the unit responsibly at the end of its working life, as it contains WEEE waste.




Vaccine Storage Rooms Installation and Operation Manual

9. LAY-OUT OF WIC 10



A - A




Mark	Qty.	Description	Changed	Date
Drawn	TTa	Order	Parts	550679-01
Drawn date	2014-12-15	Project	List	
Scale	1:25	Room		

PORKKA

PORKKA FINLAND OY
 PORKKA OY
 FIN-33101 Tampere Finland
 Pub. 020 5555 12
 Tel. +358 20 5555 12
 Fax +358 20 5555 360
 www.houffe.com

Vers. Index




Base evenness requirement $\pm 3\text{mm}/3\text{m}$ and $\pm 5\text{mm}/5\text{m}$.

This document must not be copied without our written permission and the contents thereof must not be imparted to a third party. It may be used for any unauthorized purpose.


Vaccine Storage Rooms Installation and Operation Manual

10. LAY-OUT OF WIC 20




Vaccine Storage Rooms Installation and Operation Manual

11. LAY-OUT OF WIC 30



Vaccine Storage Rooms Installation and Operation Manual

12. LAY-OUT OF WIC 40



Vaccine Storage Rooms Installation and Operation Manual

Huurre in Finland is an internationally recognised designer and manufacturer of Porkka insulated panels for cold room construction.

Huurre in Finland is an internationally renowned designer and manufacturer of insulated panels for cold room construction. As well as insulated doors, both hinged and sliding, marine fire retardant doors and specialist insulated enclosures with access hatches for food production. Clean rooms for the pharmaceutical and electronics industries is another sector where we excel.

Our factory is based in Ylöjärvi, Finland and has 14,000m² dedicated to insulated panel manufacturing and where our specialist doors and enclosures also are produced.

Huurre has over 70 years' experience in both design and production and has supplied rooms worldwide for numerous applications. Huurre will continue to invest in both research and development to ensure our customers purchase efficient, cost effective and environmentally friendly products.

Porkka is a strong brand of Huurre, well known for both supply and service in the refrigeration sector.



Huurre designs and manufactures to a quality standard ISO 9001 which is controlled and certified by Lloyd's Register Quality Assurance. We also hold the environmental certificate ISO 14001 which also has been issued by LRQA, where environment impact in production is kept to a minimum.

CE marked products. The manufacturer through continuous research and development reserves the right to change technical specification and design without notice.

Huurre reserves the right to make any changes without prior notice.



PORKKA
VALUE-ADDED REFRIGERATION

VDW CoolSystems B.V.
De Aaldor 8C
4191 PC Geldermalsen
Tel. +31 (0)345 582 546
info@vdwcoolsystems.nl
www.vdwcoolsystems.nl