

Key features

Key features

- Modular approach 1 by 1 units with corridor between
- Very easy to expand for more capacity without any disturbance on the existing freezer setup during installation of new units
- Temperature range from positive down to -80°C
- Storage- and blast freeze in one unit as well separate units
- Storage capacity up to 3,200 kg. fluid product per unit
- Extreme low deviation +/-1°C
- Full redundancy
- Low power consumption








Key features


- High personal safety, no entrance in the freezer
- Fire Rating FR60
- More than 150 alarms on each freezer unit
- Service and maintenance during full operation
- Software complies with GMP guidelines
- Factory Acceptance Test performed together with the client at our workshop in Denmark
- Very comprehensive documentation according to GDP
- Fast delivery time and limited on-site installation time
- Every installation finalizes with a Site Acceptance Test





Overview - a typical Sales drawing







Technical specifications


Technical specifications - In general

- 220 mm. FR60 rated PIR insulation panels
- FR60 rated door with 1,200 mm. door opening
- A complete closed clean room environment
- Water cooled compressor skids
- 2 pcs. 19" Siemens touch operator HMI
- Lowenco invented duct evaporator system
- Forced air ventilation method
- All components accessible from outside no maintenance people will enter the storage room



Technical specifications – Storage units


- 3 or 4 doors version
- Capacity; up to 2,400 kg. or 3,200 kg. fluid products
- Blast freeze capability can be included
- Double duct evaporator system
- Circulating the air by 3,000 m3/hour 6,000 m3/hour in assist mode
- 4 door version dimensions without tech. space:
 - L: 6,300 mm. D: 2,100 mm. H: 3,600 mm.
 - L: 5,860 mm. D: 1,660 mm. H: 3,160 mm.





Technical specifications – Blast units

- 1 or 2 door version
- Capacity; up to 800 kg. or 1,600 kg. fluid products
- Compressor technology no Liquid Nitrogen
- Full control of blast process and time
- Fast blast time e.g. 160 kg. from +8°C to -60°C in less than 4 hours and 1,600 kg. from +25°C to -60°C in less than 21 hours (custom made)
- Dimensions without tech. space:
 - L: 3,890 mm. D: 2,030 mm. H: 3,600 mm.
 - L: 3,450 mm. D: 1,590 mm. H: 3,160 mm.





Technical space







Technical space











Technical space







On-site installation







Technical details - only Best in Class components used

- Siemens S7-1500 PLC
- Siemens 19" touch HMI
- Bitzer compressors
- Danfoss VLT
- Grundfos pumps
- Fischer Kälte Klima
- ABB
- Belimo
- PR electronics temperature transmitters
- Kingspan FR60 panels



Storage systems


Trolley designs


- We design the trolley
 - to fit your products and needs bags, bottles or
 - to get the best possible utilization of capacity
- As standard two trollies goes into each door for easier handling and due to weight
- Maximum dimensions of a standard trolley: L: 1,500 mm. W: 550 mm. H: 2,140 mm.
 -but as mention, we design the trolley to fit your specific needs and wishes



Trolley designs - examples










Drawers system

- Capacity up to 800 kg. per door
- Maximum load of 100 kg. in each drawer
- Full pullout drawers
- From 8 to 16 drawers in each door
- Drawers in three different heights; 158 mm., 225 mm. and 358 mm. Inside: L: 1.405 mm. W: 535 mm.
- Height of rack is 1,700 mm.
- Drawers sizes can be mixed in each door rack
- Individual, adjustable and flexible interior design
- Location Management System integrated





Individual and full flexible interior design of drawers







Why LSSU – and not a traditional refrigeration solutions

Competitors traditional refrigeration solutions – Up-rights freezer units

- Small capacity per unit
- Large required footprint for large capacity needs
- Requires a clean room
- Air cooled compressor system
- Requires a large HVAC system
- Validation of every single freezer unit
- >85% higher power consumption
- High deviation even in such a small chamber typical >+/-5°C
- No redundancy
- No fire rating
- Not possible to get the products insured

- Limited numbers of alarms
- Limited operator information's and functions
- No service and maintenance during operation
- No Factory Acceptance Test
- Very limited documentation
- Documentation not according to GDP
- 3½ times longer pull-down time from ambient to 70°C empty
- No blast freeze capability
- Reloading and handling of all products
- Lower unit price but much higher installation cost
- Nearly ten times higher yearly Total Cost of Operation



Competitors traditional refrigeration solutions - Larger cold rooms

- Constructed on site no Factory Acceptance Test is possible
- Long construction and installation time on-site
- Safety risk when persons entering the freezer room
- Use of LN2 for blast freezing no control of the blast freezing process at all. -110°C blast into the room
- Often very limited expansion possibilities if some, it will often be very difficult
- High deviation in the large room typical >+/-5°C





