












EXTERIOR

A modern design in combination with an extremely strong box frame construction. Finished with a scratch proof epoxy coating, colour grey-white RAL 9010. Standard mounted on adjustable legs.

INTERIOR

The interior is made of high quality, easy to clean stainless steel (DIN 1.4301). The use of a tangential fan in the return air maintains the highest temperature uniformity possible. A lockable entry port, water drain and height adjustable stainless steel

Fig. 1 Aircirculation



CONDITIONING

A powerful, but energy saving, climate system for heating and cooling guarantees a high temperature stability and fast temperature recovery. The tangential fan in the return air guarantees the temperature uniformity. The use of this type of fan helps ensure that the same conditions are reached at every point within the chamber removing cold spots and hence condensation within the entire chamber

FLEXIBILITY


All incubators can be delivered with or without cooling. All models can be installed side by side to form an incubator wall. Practical division for each application within the chamber is also possible, e.g by height adjustable platforms or a drawer system mounted on telescopic arms.

EFFECTIVITY

Sustainability, environmental friendly and lower energy usage were the keyword issues in the development of the incubator range. In other words as low a carbon footprint as possible. This for example has been achieved by using durable material sourced from European suppliers only and efficient cooling using a bypass



Fig. 2 Configuration: linear calibration




system, driven by the latest energy saving compressors. The slim line design optimizes capacity within the chamber while maintaining a small footprint on the laboratory floor.

TEMPERATURE CONTROLLER


The controller displays temperature and time digitally and can be programmed in 8 Temperature Setpoints (TS-1 upto TS-8) in association with 8 Time-units (T-1 upto T-8). The temperature is linear and hence ensures the tightest possible control of temperature over the entire temperature range (Fig. 3) By installing additional software it is also possible to link multiple incubators together for remote monitoring and control. External temperature alarms are routed via a potential-free contact or a SMS alert system over the internet.

SAFETY

Various safety features have been built-in for the protection of the product within the incubator. Firstly there is an alarm activated from within the controller if the band width is exceeded by +/- 5°C from set. Secondly there are backup safety thermostats for both high and low temperatures. All alarms are acoustic and optical.



Program example





AUTO DEFROST

In durability tests longer than 24 hours and below +4°C, an auto defrost has to be used to avoid ice forming on the evaporator. The frequency, temperature and length of time for the auto defrost system is adjustable. The settings for the defrost cycle(s) are dependent on the required set temperature for the chamber and the usage of the incubator.

ANTI-CONDENSATION SYSTEM

Desiccation occurs when humidity is drawn, for example, from the media in a Petri dish due to a "cold spot" where it condenses out of the atmosphere. The unique Snijders anti condensation system eliminates the potential for "cold spots" inside the working area and hence avoids desiccation. This system consigns the problem of desiccation to history.

(COOLED) INCUBATORS AND WALLS

TECHNICAL SPECIFICATIONS

CRIPTION : volume ss volume erial internal	INCUBATOR BR1 494 liter	INCUBATOR BR2	COOLED INCUBATOR EB1-N / EB1-D	INCUBATOR
volume ss volume erial internal	BR1	BR2		
volume ss volume erial internal			EB1-N / EB1-D	EDO NIZEDO D
ss volume erial internal	494 liter	(0.) 074 !!!		EB2-N / EB2-D
erial internal		(2x) 231 liter	494 liter	(2x) 231 liter
	576 liter	(2x) 269 liter	576 liter	(2x) 269 liter
	Stainless steel	Stainless steel	Stainless steel	Stainless steel
erial external	Coated steel	Coated steel	Coated steel	Coated steel
lation	Polyurethane	Polyurethane	Polyurethane	Polyurethane
ernal dimensions, (wxdxh)	820x855x1960 mm	820x855x1960 mm	820x855x1960 mm	820x855x1960 mm
rnal dimensions (wxdxh)	635x600x1295 mm	(2x) 635x600x605 mm	635x600x1295 mm	(2x) 635x600x605 mm
form dimensions, (bxd)	625x580 mm	625x580 mm	625x580 mm	625x580 mm
mal transportation (wxdxh)	820x775x1945 mm	820x775x1945 mm	820x775x1945 mm	820x775x1945 mm
forms standard /	5/29	(2x) 2/10	5/29	(2x) 2/10
. quantity				
ance between platforms	40 mm	40 mm	40 mm	40 mm
y port Ø	45 mm	45 mm	45 mm	45 mm
er within controller	Digital	Digital	Digital	Digital
troller	Electronic PID	Electronic PID	Electronic PID	Electronic PID
perature sensor	PT100	PT100	PT100	PT100
perature safety	Adjustable by high	Adjustable by high	Adjustable by high	Adjustable by high
	temperature	temperature	temperature	temperature
	thermostat with	thermostat with	thermostat with	thermostat with
	separate PT100 sensors	separate PT100 sensors	separate PT100 sensors	separate PT100 sensors
	No	No	Optional	Optional
troller perature sensor	Electronic PID PT100 Adjustable by high temperature thermostat with separate PT100 sensors	Electronic PID PT100 Adjustable by high temperature thermostat with separate PT100 sensors	Electronic PID PT100 Adjustable by high temperature thermostat with separate PT100 sensors	Electronic PID PT100 Adjustable by temperature thermostat wi separate PT10

SPECIFICATIONS (DIN 12880)						
Temperature range, in °C	(Tamb. +8°C) tot +60	(Tamb.+8°C) tot +60	+4 tot +60 / -10 tot +60	+4 tot +60 / -10 tot +60		
Temperature setting, in °C	0,1	0,1	0,1	0,1		
Temperature fluctuation, in °C	≤ 0,2 (at 1 spot)	≤ 0,2 (at 1 spot)	≤ 0,2 (at 1 spot)	≤ 0,2 (at 1 spot)		
Temperature variation, in °C	\leq 0,5 (in chamber)	\leq 0,5 (in chamber)	\leq 0,5 (in chamber)	\leq 0,5 (in chamber)		
Airflow m/s	0,2	0,2	0,2	0,2		

Left or right	Left or right	Left or right	Left or right
(at the right= ex works)	(at the right= ex works)	(at the right= ex works)	(at the right= ex works)
Grip incl. cylinder lock	Grip incl. cylinder lock	Grip incl. cylinder lock	Grip incl. cylinder lock
30 x 40 mm, adjustable,	30 x 40 mm, adjustable,	30 x 40 mm, adjustable,	30 x 40 mm, adjustable,
height 150 mm	height 150 mm	height 150 mm	height 150 mm
220-240V; 50 Hz.	220-240V; 50 Hz.	220-240V; 50 Hz.	220-240V; 50 Hz.
196	209	224	237
	(at the right= ex works) Grip incl. cylinder lock 30 x 40 mm, adjustable, height 150 mm 220-240V; 50 Hz.	(at the right= ex works) (at the right= ex works) Grip incl. cylinder lock Grip incl. cylinder lock Grip incl. cylinder lock 30 x 40 mm, adjustable, height 150 mm height 150 mm 220-240V; 50 Hz. 220-240V; 50 Hz.	(at the right= ex works) (at the right= ex works) (at the right= ex works) Grip incl. cylinder lock Height 150 mm, adjustable, 30 x 40 mm, adjustable, height 150 mm height 150 mm 220-240V; 50 Hz. 220-240V; 50 Hz.

^{*}Specifications subject to change

SNIJDERS LABS: EXPERIENCED INNOVATORS

SNIJDERS LABS forms part of the Snijders Group, which actively delivers equipment and products for scientific research & development as well as internal transport systems and examination couches for the health care sector under SNIJDERS CARE. All design, manufacture and testing is held in house to assure high quality production and investment in new technologies for the production of all Snijders products. The total control of the manufacturing line means that Snijders can offer total quality, in-depth knowledge and detailed assistance to all of their clients.

Tweede druk 2014















SERVICE AND WARRANTY

Contact your local distributor who will guarantee quality and service (if necessary check our website for distributor details).

VISIT OUR WEBSITE WWW.SNIJDERSLABS.COM

There you'll find all the latest information about:

- + ULT freezers (-86°C) with datasheets of any type, racking systems, boxes and other accessories
- a variety of climate cabinets for plants, seed germination, fungi, snails and insects research with temperature-, light- and humidity control
- + (cooled) incubators and incubator walls, designed for general microbiological research of among others food, water and medical laboratories.

