



# Operating Instructions

Ultra-Low Temperature Freezer

## MDF-1156

## MDF-1156ATN Series



MDF-1156

Please read the operating instructions carefully before using this product, and keep the operating instructions for future use.

See page 30 for all model numbers.



# CONTENTS

|                                |       |
|--------------------------------|-------|
| INTRODUCTION                   | P. 2  |
| PRECAUTIONS FOR SAFE OPERATION | P. 3  |
| ENVIRONMENTAL CONDITIONS       | P. 7  |
| FREEZER COMPONENTS             | P. 8  |
| Control panel                  | P. 10 |
| INSTALLATION SITE              | P. 11 |
| INSTALLATION                   | P. 12 |
| START-UP OF UNIT               | P. 13 |
| CHAMBER TEMPERATURE SETTING    | P. 14 |
| REMOTE ALARM TERMINAL          | P. 14 |
| ALARM TEMPERATURE SETTING      | P. 15 |
| SETTING OF ALARM RESUME TIME   | P. 16 |
| ALARMS & SAFETY FUNCTIONS      | P. 17 |
| ROUTINE MAINTENANCE            | P. 18 |
| Cleaning of cabinet            | P. 18 |
| Cleaning of condenser filter   | P. 18 |
| Defrosting of inside wall      | P. 19 |
| TROUBLESHOOTING                | P. 20 |
| DISPOSAL OF UNIT               | P. 21 |
| Recycle of battery             | P. 21 |
| Decontamination of unit        | P. 21 |
| DISPOSAL OF BATTERY            | P. 26 |
| TEMPERATURE RECORDER           | P. 27 |
| BACKUP COOLING KIT             | P. 28 |
| SPECIFICATIONS                 | P. 29 |
| PERFORMANCE                    | P. 30 |
| SAFETY CHECK SHEET             | P. 31 |

# INTRODUCTION

- Read the operating instructions carefully before using the product and follow the instructions for safe operation.
- PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.
- Keep the operating instructions in a suitable place so that they can be referred to as necessary.
- The operating instructions are subject to change without notice for improvement of performance or function.
- Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.
- No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

## **IMPORTANT NOTICE**

PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for any loss or damage to the contents of the product.

# PRECAUTIONS FOR SAFE OPERATION

**It is imperative that the user complies with the operating instructions as it contains important safety advice.**

Items and procedures are described so that you can use this unit correctly and safely. If the precautions advised are followed, this will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:




## **WARNING**

Failure to observe **WARNING** signs could result in a hazard to personnel possibly resulting in serious injury or death.

## **CAUTION**

Failure to observe **CAUTION** signs could result in injury to personnel and damage to the unit and associated property.

Symbol shows;

-  This symbol means caution.
-  This symbol means an action is prohibited.
-  This symbol means an instruction must be followed.

Be sure to keep the operating instructions in a place accessible to users of this unit.

< Label on the unit >



This mark is labeled on the cover in which the electrical components of high voltage are enclosed to prevent the electric shock.















The cover should be removed by a qualified engineer or a service personnel only.

## **WARNING**

As with any equipment that uses CO<sub>2</sub> or N<sub>2</sub> gas, there is a likelihood of oxygen depletion in the vicinity of the equipment. It is important that you assess the work site to ensure there is suitable and sufficient ventilation. If restricted ventilation is suspected, then other methods of ensuring a safe environment must be considered. These may include atmosphere monitoring and warning devices.


# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Do not use the unit outdoors.** Current leakage or electric shock may result if the unit is exposed to rain water.
-  **Only qualified engineers or service personnel should install the unit.** The installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit on a sturdy floor and take an adequate precaution to prevent the unit from turning over.** If the floor is not strong enough or the installation site is not adequate, this may result in injury from the unit falling or tipping over.
-  **Never install the unit in a humid place or a place where it is likely to be splashed by water.** Deterioration of the insulation may result which could cause current leakage or electric shock.
-  **Never install the unit in a flammable or volatile location.** This may cause explosion or fire.
-  **Never install the unit where acid or corrosive gases are present** as current leakage or electric shock may result due to corrosion.
-  **Always ground (earth) the unit to prevent electric shock.** If the power supply outlet is not grounded, it will be necessary to install a ground by qualified engineers.
-  **Never ground the unit through a gas pipe, water main, telephone line or lightning rod.** Such grounding may cause electric shock in the case of an incomplete circuit.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **Never store volatile or flammable substances** in this unit if the container cannot be sealed. These may cause explosion or fire.
-  **Do not insert metal objects such as a pin or a wire into any vent, gap or any outlet on the unit.** This may cause electric shock or injury by accidental contact with moving parts.
-  **Use this unit in safe area when treating the poison, harmful or radiate articles.** Improper use may cause bad effect on your health or environment.
-  **Turn off the power switch (if provided) and disconnect the power supply to the unit prior to any repair or maintenance** of the unit in order to prevent electric shock or injury.
-  **Do not touch any electrical parts (such as power supply plug) or operate switches with a wet hand.** This may cause electric shock.








# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Ensure you do not inhale or consume medication or aerosols** from around the unit at the time of maintenance. These may be harmful to your health.
-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never put cylinders with liquid on the unit** as this may cause electric shock or short circuit when the liquid is spilled.
-  **Never bind, process, or step on the power supply cord, or never damage or break the power supply plug.** A broken supply cord or plug may cause fire or electric shock.
-  **Do not use the supply cord if its plug is loose.** Such supply cord may cause fire or electric shock.
-  **Never disassemble, repair, or modify the unit yourself.** Any such work carried out by an unauthorized person may result in fire, or electric shock or injury due to a malfunction.
-  **Disconnect the power supply plug if there is something wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  **When removing the plug from the power supply outlet, grip the power supply plug, not the cord.** Pulling the cord may result in electric shock or fire by short circuit.
-  **Disconnect the power supply plug** before moving the unit. Take care not to damage the power cord. A damaged cord may cause electric shock or fire.
-  **Disconnect the power plug when the unit is not used for long periods.** Keeping the connection may cause electric shock, current leakage, or fire due to the deterioration of insulation.
-  **If the unit is to be stored unused in an unsupervised area for an extended period, ensure that children do not have access and that doors cannot be closed completely with a key.**
-  **The disposal of the unit should be accomplished by appropriate personnel.** Remove doors to prevent accidents such as suffocation.
-  **Do not put the packing plastic bag within reach of children** as suffocation may result.
-  **Fix the shelf stoppers and shelves securely.** Incomplete installation may cause injury or damage.
-  **Do not place this unit in a location where it is difficult to disconnect the power supply plug.** Failure to disconnect the power supply plug may cause fire in the event of a problem or malfunction.

# PRECAUTIONS FOR SAFE OPERATION

## CAUTION

-  This unit must be plug into a dedicated circuit protected by branch circuit breaker.
-  Use a dedicated power source as indicated on the rating label attached to the unit. A multiple-tap may cause fire resulting from abnormal heating.
-  **Connect the power supply plug to the power source firmly after removing the dust on the plug.** A dusty plug or improper insertion may cause a heat or ignition.
-  **Never store corrosive substances such as acid or alkali** in this unit if the container cannot be sealed. These may cause corrosion of inner components or electric parts.
-  **Check the setting when starting up of operation after power failure or turning off of power switch.** The stored items may be damaged due to the change of setting.
-  **Be careful not to tip over the unit** during movement to prevent damage or injury.
-  **Prepare a safety check sheet** when you request any repair or maintenance for the safety of service personnel.




# ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):

- Indoor use;
- Altitude up to 2000 m;
- Ambient temperature 5°C to 40°C;
- Maximum relative humidity 80% for temperature up to 31°C decreasing linearly to 50% relative humidity at 40°C;
- Mains supply voltage fluctuations up to  $\pm 10\%$  of the nominal voltage;
- Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;
- Temporary OVERVOLTAGES occurring on the mains supply;
- Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases)

# FREEZER COMPONENTS



MDF-1156ATN

# FREEZER COMPONENTS

**1. Lock**

**2. Door:** Hinged type. The door can be opened in any angle on the way to full open.

**3. Magnetic door gasket:** Seals the door and prevents leakage of cold air.

**4. Inner lid:** Serves as a means of reducing cold air leakage when the door is open.

**5. Access port:** Serves a means of leading the measuring cable from the freezer chamber to the outside.

**6. Caster:** 6 casters are provided. They make the moving of the unit easier.

**7. Leveling foot:** Serves to adjust the height and to settle the frame.

**8. Grille:** Acts as an inlet for air to cool the motor. Be careful not to block this. By pulling down this grille, you can clean a clogged condenser filter.

**9. Exhaust air vent:** Be careful not to block this.

**10. Control panel cover lock:** To lock the control panel cover to avoid the setting by accidental contact.

**11. Control panel cover**

**12. Control panel:** Refer to page 10.

**13. Digital temperature indicator:** This indicator shows the present temperature or setting temperature.


**14. Temperature recorder (ATN type)**

**15. Remote alarm terminal:** It is located inside of rear cover. Refer to page 14 "Remote alarm terminal".

**16. Backup cooling kit joint (ATN type):** It is positioned at rear of the unit. Serves to connect with the pipelines from the cylinder (liquid N<sub>2</sub>) at the top right of the rear frame. Refer to page 28 "Backup cooling kit"

# FREEZER COMPONENTS

## Control panel



<MDF-1156ATN>

- 1. Power supply switch (POWER):** Power all functions except remote alarm and backup cooling kit for ATN type.
- 2. Backup switch (BACK UP) (ATN type):** Switch ON for operation of the backup cooling kit and switch OFF for stopping. Handle this switch according to page 28 "Backup cooling kit".
- 3. Backup test switch (TEST) (ATN type):** Examine the functions of the backup cooling kit. This switch allows liquid N<sub>2</sub> to spout under any circumstances. Handle it according to page 28 "Backup cooling kit".
- 4. Remote alarm switch (REMOTE ALARM):** This switch is for remote alarm. In case of operating remote alarm, turn the switch to ON.
- 5. Alarm buzzer stop key (BUZZER):** To stop the alarm from sounding, press this key. Should a further abnormality occur, the buzzer will sound automatically.
- 6. Alarm test key (ALARM):** Check that the alarm lamp and the buzzer are functional when the freezer is operating well.
- 7. Temperature setting key (PV/SV):** This key has two functions.  
 PV; present freezer compartment temperature is displayed. (+50 °C to -170°C)  
 SV; setting temperature is displayed. (-130°C to -152°C)  
 When the desired temperature is set, this key should be pressed (SV side). When no key is pressed more than 90 seconds, displayed number won't be memorized, setting will be over and digital temperature indicator will display the present value.
- 8. Digit shift key (▶▶):** The digit of the figure displayed for temperature adjustment can be shifted using this key.
- 9. Numerical value shift key (▲):** The figure that is displayed digitally can be changed by pressing this key.
- 10. Enter key (ENT):** By pressing this key in set mode, the setting is memorized and set mode is completed.
- 11. Alarm lamp (ALARM):** This lamp will flash when the unit is in alarm condition.
- 12. Filter check lamp (FILTER CHECK):** This lamp blinks and the alarm buzzer sounds when the condenser filter is clogged. Clean the condenser filter according to page 18 "Cleaning of condenser filter".
- 13. Backup lamp:** This lamp will light when the backup switch is ON. (This doesn't show the case that backup cooling kit is activated.)
- 14. Digital temperature indicator (TEMPERATURE °C):** This indicator shows the present temperature or setting temperature (flash).
- 15. Temperature recorder (ATN type):** This records the internal temperature. Refer to page 27 "Temperature recorder".
- 16. Battery switch:** It is the switch of the battery for the power failure alarm. Turn it on usually.

# INSTALLATION SITE

To operate this unit properly and to obtain maximum performance, install the unit in a location with the following conditions:

■ **A location not subjected to direct sunlight**

Do not install the unit under direct sunlight. Installation in a location subjected to direct sunlight cannot obtain the intended performance.

■ **A location with adequate ventilation**

Leave at least 10 cm around the unit for ventilation. Poor ventilation will result in a reduction of the performance and consequently the failure.

■ **A location away from heat generating sources**

Avoid installing the unit near heat-emitting appliances such as a heater or a boiler etc. Heat can decrease the intended performance of the unit.

■ **A location with little temperature change**

Install the unit under stable ambient temperature. The allowable ambient temperature is between -5 and +35°C.

■ **A location with a sturdy and level floor**

Always install the unit on a sturdy and level floor. The uneven floor or tilted installation may cause failure or injury. Install the unit in stable condition to avoid the vibration or noise. Unstable condition may cause vibration or noise.

 **WARNING**

**Install the unit on a sturdy floor.** If the floor is not strong enough or the installation site is not adequate, this may result in injury from the unit falling or tipping over.

**Select a level and sturdy floor for installation.** This precaution will prevent the unit from tipping. Improper installation may result in water spillage or injury from the unit tipping over.

■ **A location not prone to high humidity**

Install the unit in the ambient of 80% R.H. or less humidity. Installation under high humidity may cause current leakage or electric shock.

 **WARNING**

**Do not use the unit outdoors.** Current leakage or electric shock may result if the unit is exposed to rain water.

**Never install the unit in a humid place or a place where it is likely to be splashed by water.** Deterioration of the insulation may result which could cause current leakage or electric shock.

■ **A location without flammable or corrosive gas**

Never install the unit in a flammable or volatile location. This may cause explosion or fire or may result in the current leakage or electric shock by the corrosion of the electrical components.

■ **A location without the possibility of anything fall**

Avoid installing the unit in the location where anything can fall down onto the unit. This may cause the breakdown or failure of the unit.

# INSTALLATION

## 1. Removing the packaging materials and tapes


Remove all transportation packaging materials and tapes. Open the doors and ventilate the unit. If the outside panels are dirty, clean them with a diluted neutral dishwashing detergent. (Undiluted detergent can damage the plastic components. For the dilution, refer to the instruction of the detergent.) After the cleaning with the diluted detergent, always wipe it off with a wet cloth. Then wipe off the panels with a dry cloth.

### Note:

Remove the cable tie banding the power supply cord. Prolonged banding may cause the corrosion of the cord coating.

## 2. Adjusting the leveling foot

Extend the leveling feet by rotating them counterclockwise to contact them to the floor. Ensure the unit is level.



## 3. Fixing the unit

Two fixtures are attached to the rear of the frame. Fix the frame to the wall with these fixtures and rope or chain.

## 4. Ground (earth)

The ground (earth) is for preventing the electric shock in the case of the electrical insulation is somehow degraded. Always ground the unit at the time of installation.

### **⚠ WARNING**

**Use a power supply outlet with ground (earth) to prevent electric shock.** If the power supply outlet is not grounded, it is necessary to install a ground by qualified engineers.

**Never ground the unit through a gas pipe, water main, telephone line or lightning rod.** Such grounding may cause electric shock in the case of an incomplete circuit.

# START-UP OF UNIT

---

Follow the procedures for the initial and consequent operations of the unit.







1. Make sure that all the switches on the control panel, such as the power supply switch, battery switch, backup switch (ATN type) and remote alarm switch are off.
2. Connect the power cord to the dedicated outlet having appropriate rating with the chamber empty, and turn on the power supply switch on the freezer.
3. Turn on the battery switch.
4. The audible alarm may be activated. In this case, press the alarm buzzer stop key (BUZZER) to silence the alarm.
5. Set the desired chamber temperature. See page 14 for the chamber temperature setting.
6. Check that the chamber temperature reaches the desired temperature.
7. Turn on the backup switch (ATN type) and remote alarm switch.
8. Make sure that the alarm lamp blinks and the buzzer sounds by pressing the alarm test key (ALARM).
9. After confirming the above, you can put articles into the freezer chamber in a small batch to prevent the temperature rise.

# CHAMBER TEMPERATURE SETTING

Table 1 shows the basic procedure for setting the chamber temperature. Perform key operations in the sequence indicated in the table. The example in the table is based on the assumption that the desired temperature is -145°C.

**Note:** The unit is set at the factory that the chamber temperature -150°C.

**Table 1. Basic operation sequence (Example: Internal temperature -145°C)**

|   | Description of operation                                      | Key operated | Indication after operation                                                                                                                                         |
|---|---------------------------------------------------------------|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Turn on the power switch.                                     | ----         | The third digit blinks.                                                         |
| 2 | Set the second digit with the digit shift key.                | ▶▶           | The second digit blinks.                                                        |
| 3 | Set the second digit to 4 with the numerical value shift key. | ▲▲           | When pressed, the figure of settable digit changes.                             |
| 4 | Set the first digit with the digit shift key.                 | ▶▶           | When pressed, the settable digit is shifted.                                    |
| 5 | Set the first digit to 5 with the numerical value shift key.  | ▲▲           | When pressed, the figure of settable digit changes.                             |
| 6 | Press enter key.                                              | ENT          | Set temperature is memorized and the current chamber temperature is displayed.  |

**Note:**


- The temperature set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.
- Although the value of the chamber temperature setting can range from -125°C to -155°C, the guaranteed temperature with no load is -152°C when the ambient temperature is 30°C.

# REMOTE ALARM TERMINAL

The terminal of the remote alarm is installed at the back of the unit.

The signal is contact output. Contact capacity is 2 A (DC 30 V)/0.1 mA (DC 1 V).

- a) output: normal open, connect with N.O. and COM.
- b) output: normal close, connect with N.C. and COM.





# ALARM TEMPERATURE SETTING

This unit is provided with the high temperature alarm. The setting of high temperature alarm is 10°C or 15°C (The setting value is 2 kinds.) higher than the setting of chamber temperature. The procedure in table 2 shows the sequence to set the high temperature alarm at 15°C higher than the setting of chamber temperature. Set the figure to 000 at step 5 when setting the high temperature alarm at 10°C higher than the setting of the set temperature.

**Note:**

The alarm temperature is set at the factory 10°C higher than the setting of chamber temperature.

**Table 2. High temperature alarm setting (15°C higher than the setting of chamber temperature)**

|   | Description of operation                                    | Key operated | Indication after operation                                                       |
|---|-------------------------------------------------------------|--------------|----------------------------------------------------------------------------------|
| 1 |                                                             | ----         | Set temperature is displayed.                                                    |
| 2 | Press the numerical value shift key for about 5 seconds.    |              | The first digit blinks.                                                          |
| 3 | Press numerical value shift key and scroll the figure to 1. |              | The first digit blinks.                                                          |
| 4 | Press enter key.                                            | ENT          | The first digit blinks.                                                          |
| 5 | Scroll the figure to 1 by using numerical value shift key.  |              | When pressed, the figure of settable digit changes.                              |
| 6 | Press enter key.                                            | ENT          | Alarm temperature is memorized and the current chamber temperature is displayed. |

**Note:** The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

# SETTING OF ALARM RESUME TIME

The alarm buzzer is stopped by pressing alarm buzzer stop key (BUZZER) on the control panel during alarm condition (Initial setting 130: The remote alarm is not stopped).

The buzzer will be activated again after certain suspension if the alarm condition is continued. The suspension time can be set by following the procedure shown in the Table 3 below.

The example in the table is based on the assumption that the desired duration is 20 minutes.

**Note:** The duration is set to 130 at the factory.

**Table 3. Setting procedure (change to 20 minutes from no alarm resume time)**

|   | Description of operation                                                         | Key operated | Indication after operation                                                       |
|---|----------------------------------------------------------------------------------|--------------|----------------------------------------------------------------------------------|
| 1 |                                                                                  | ----         | The current chamber temperature is displayed.                                    |
| 2 | Press numerical value shift key for about 5 seconds.                             |              | The first digit blinks.                                                          |
| 3 | Set the figure to F25 with the digit shift key and numerical value shift key.    |              | The settable digit is shifted.                                                   |
|   |                                                                                  |              | When pressed, the figure of settable digit changes.                              |
| 4 | Press enter key.                                                                 | ENT          | The current setting is displayed and the first digit blinks.                     |
| 5 | Scroll the figure to 120 with the digit shift key and numerical value shift key. |              | The settable digit is shifted.                                                   |
|   |                                                                                  |              | When pressed, the figure of settable digit changes.                              |
| 6 | Press enter key.                                                                 | ENT          | Alarm resume time is memorized and the current chamber temperature is displayed. |

- The settable alarm resume time is between 10 and 60 minutes with 1 minute increment. (The setting is 110, 111, ----- 159, 160). The buzzer would not resume if the resume time is set in 100.
- The buzzer is ON when the resume time is changed in F25 while the resume time is counted (the buzzer is silenced by alarm buzzer stop key (BUZZER). The count of resume time is stopped.
- The setting cannot be changed during power failure.
- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing enter key (ENT) is not memorized.
- It is selectable whether the buzzer is in conjunction with the remote alarm or not when the buzzer is silenced by alarm buzzer stop key (BUZZER). (The remote alarm is not in conjunction with the buzzer by factory setting.)
- The remote alarm is activated in conjunction with the buzzer if the third digit "1" is set to "0" in the step 5 above

# ALARMS & SAFETY FUNCTIONS

This unit has the alarm and safety functions shown in table below, and also self diagnostic function.

**Table 4. Alarms and safety functions**

| Alarm & safety              | Situation                                                                  | Indication                                                          | Buzzer                               | Safety operation                                      |
|-----------------------------|----------------------------------------------------------------------------|---------------------------------------------------------------------|--------------------------------------|-------------------------------------------------------|
| High temperature alarm      | When the chamber temp. is approx. 15°C (or 10°C) higher than the set temp. | Alarm lamp blinks.                                                  | Intermittent tone with 12 min. delay | Remote alarm with 12 min. delay                       |
| Power failure alarm         | At power failure<br>When power cord is disconnected.                       | Alarm lamp blinks.                                                  | Intermittent tone                    | Remote alarm                                          |
| Filter check                | When the condenser filter is clogged.                                      | Filter check lamp blinks.                                           | Intermittent tone                    | -----                                                 |
| Auto-return                 | When there is no key pressing in each setting mode for 90 seconds.         | Chamber temp. is displayed.                                         | -----                                | Finishing of each setting mode                        |
| Sensor abnormality          | If the thermal sensor is disconnected.                                     | Alarm lamp blinks. E01 and 50°C is displayed alternately.           | Intermittent tone                    | Remote alarm<br>Continuous operation                  |
|                             | If the thermal sensor is short-circuited.                                  | Alarm lamp blinks. E02 and -170°C are displayed alternately.        |                                      |                                                       |
|                             | If the cascade sensor is disconnected.                                     | Alarm lamp blinks. E03 and chamber temp. are displayed alternately. | Intermittent tone                    | Remote alarm                                          |
|                             | If the cascade sensor is short-circuited.                                  | Alarm lamp blinks. E04 and chamber temp. are displayed alternately. |                                      |                                                       |
|                             | If the filter sensor is disconnected.                                      | Alarm lamp blinks. E05 and chamber temp. are displayed alternately. | Intermittent tone                    | Remote alarm                                          |
|                             | If the filter sensor is short-circuited.                                   | Alarm lamp blinks. E06 and chamber temp. are displayed alternately. |                                      |                                                       |
|                             | If the ambient temperature sensor is disconnected.                         | ALARM lamp blinks. E07 and chamber temp. are displayed alternately. | Intermittent tone                    | Remote alarm                                          |
|                             | If the ambient temperature sensor is short-circuited.                      | ALARM lamp blinks. E08 and chamber temp. are displayed alternately. |                                      |                                                       |
| Condenser temp. abnormality | In the event of failure of fan motor for cooling the compressor            | E10 and chamber temp. are displayed alternately.                    | Intermittent tone                    | Remote alarm.<br>Compressor of high stage side stops. |

**Note:**

The filter check lamp sometimes blinks at the start-up under high ambient temperature. The lamp is off when the chamber temperature is getting lower.

**Operation after power failure**

The set value (setting of chamber temperature and alarm temperature) is memorized by nonvolatile memory. Accordingly, the freezer resumes the operation with setting before power failure.

This unit has the alarms and safety functions shown in Table 4, and also self diagnostic functions.

# ROUTINE MAINTENANCE

## **⚠ WARNING**

**Always disconnect the power supply to the unit prior to any repair or maintenance** of the unit in order to prevent electric shock or injury.

**Ensure you do not inhale or consume medication or aerosols** from around the unit at the time of maintenance. These may be harmful to your health.


## **Cleaning of cabinet**

- Clean the unit once a month. Regular cleaning keeps the unit looking new.
- Use a dry cloth to wipe off small amounts of dirt on the outside and inside of the unit and all accessories. If the outside panels are dirty, clean them with a diluted neutral dishwashing detergent. (Undiluted detergent can damage the plastic components. For the dilution, refer to the instruction of the detergent.) After the cleaning with the diluted detergent, always wipe it off with a wet cloth. Then wipe off the cabinet or accessories with a dry cloth.
- Never pour water onto or into the unit. Doing so can damage the electric insulation and cause failure.
- The compressor and other mechanical parts are completely sealed. This unit requires absolutely no lubrication.
- Check the backup cooling kit by pressing backup test switch once a month if it is installed.
- Remove the frost or ice on the chamber wall and clean the condenser filter once a month.

## **Cleaning of condenser filter**

This unit is provided with the filter check lamp. This lamp blinks and the alarm buzzer sounds when the condenser filter is clogged. Clean the condenser filter according to the following procedure. As a clogged condenser filter may cause poor cooling and compressor trouble, clean it once a month.

1. The condenser filter fitted to the grille can be removed as shown in the figure.
2. Pull out the grille.
3. Remove the condenser filter installed inside the grille and wash with water and dry it naturally.
4. Replace the clean condenser filter as before.
5. Replace the grille.



## **⚠ WARNING**

**Do not touch the condenser directly** when the condenser filter is removed for cleaning. This may cause injury by hot surface.

# ROUTINE MAINTENANCE

## Defrosting of inside wall

Defrost the inside wall of the freezer as follows:

### **Normal defrosting**

Remove the frost by the enclosed scraper.

### **Thorough defrosting**

1. Take out and transfer all the contents to another freezer or container which contains liquid CO<sub>2</sub>, liquid N<sub>2</sub>, or dry ice. Switch off the remote alarm, battery switch and backup cooling kit. Switch off the power supply switch.
2. Open the door and remove the inner lid. Leave the freezer as it is. The water remaining in the freezer compartment should be wiped up.
3. After cleaning is completed, restart the operation according to the procedure on page 13. Put back the articles into the sufficiently cooled freezer compartment.

# TROUBLESHOOTING

If the unit malfunctions, check out the following before calling for service.

| Malfunction                                                                            | Check/Remedy                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| If nothing operates even when switched on                                              | <ul style="list-style-type: none"> <li>■ The unit is not connected to the power supply.</li> <li>■ There is a power failure.</li> <li>■ The fuse is blown or the circuit breaker is activated.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                   |
| No refrigeration                                                                       | <ul style="list-style-type: none"> <li>■ The voltage too low. (In this case, call an electrician.)</li> <li>■ The breaker free.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| The cooling is poor                                                                    | <ul style="list-style-type: none"> <li>■ The environmental temperature is too high.</li> <li>■ The door is not shut tightly.</li> <li>■ The inner lid is not installed correctly.</li> <li>■ The set temperature is not set properly.</li> <li>■ The grille is blocked out.</li> <li>■ The condenser filter is clogged.</li> <li>■ The freezer is in the direct sunlight.</li> <li>■ There is any heating source near the freezer.</li> <li>■ A rubber cap and insulation for the access port are not set correctly.</li> <li>■ You put too many unfrozen articles into the freezer compartment.</li> </ul> |
| Alarm test key cannot actuate the alarm                                                | <ul style="list-style-type: none"> <li>■ The battery has been discharged entirely. In this case, operate the freezer for about 3 hours and switch on the alarm test key again. Entirely discharged cell requires about 2 days' operation.</li> <li>■ When only the buzzer or only the alarm is actuated by the alarm test key, the unactuated part is out of order, and must be replaced.</li> </ul>                                                                                                                                                                                                        |
| Recorder for ATN type does not operate normally                                        | <ul style="list-style-type: none"> <li>■ If the chart was not fed, examine the dry cell. The cell for the recorder lasts about 1 year.</li> <li>■ If the temperature was not recorded, examine the ink pen. The pen should be our specified one. Is the recorder chart set correctly?</li> </ul>                                                                                                                                                                                                                                                                                                            |
| Noise                                                                                  | <ul style="list-style-type: none"> <li>■ The freezer is not stable.</li> <li>■ The freezer is not leveled with the leveling feet.</li> <li>■ There is anything touching the frame.</li> <li>■ It is the first operation after shutdown.</li> <li>■ High temperature loads this freezer, and sometimes causes a noise. Accordingly when the temperature lowered, the noise ceases.</li> </ul>                                                                                                                                                                                                                |
| Backup test switch does not operate normally (when the backup cooling kit is provided) | <ul style="list-style-type: none"> <li>■ The liquid N<sub>2</sub> tank is not filled up.</li> <li>■ The valve of the tank is not opened.</li> <li>■ The pressure of liquid N<sub>2</sub> is 0.5 kg/cm<sup>2</sup>G. Inquire at liquid N<sub>2</sub> suppliers about its installation, removal, adjustment, and examination.</li> </ul>                                                                                                                                                                                                                                                                      |

**Note:**

If the malfunction is not eliminated after checking the above items, or the malfunction is not shown in the above table, contact our sales representative or agent.

# DISPOSAL OF UNIT

## WARNING

If the unit is to be stored unused in an unsupervised area for an extended period **ensure that children do not have access and doors cannot be closed completely.**

**The disposal of the unit should be accomplished by appropriate personnel. Always remove doors** to prevent accidents such as suffocation.

## Recycle of battery



Ni-MH

- Label indication is obliged to comply with Japanese battery regulation.



Pb

- Label indication is obliged to comply with Japanese battery regulation.



- Label indication is obliged to comply with Taiwanese battery regulation.

## Decontamination of unit

Before disposing a ultra low temperature freezer with biohazardous danger, decontaminate the ultra low temperature freezer to the extent possible by the user.

# DISPOSAL OF UNIT

(English)

## Disposal of Old Equipment and Batteries

Only for European Union and countries with recycling systems



These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.



By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.



Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

### Note for the battery symbol (bottom symbol):

This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

**Pb**

(German)

## Entsorgung von Altgeräten und Batterien

Nur für die Europäische Union und Länder mit Recyclingsystemen



Dieses Symbol, auf den Produkten, der Verpackung und/oder den Begleitdokumenten, bedeutet, dass gebrauchte elektrische und elektronische Produkte sowie Batterien nicht in den allgemeinen Hausmüll gegeben werden dürfen.

Bitte führen Sie alte Produkte und verbrauchte Batterien zur Behandlung, Aufarbeitung bzw. zum Recycling gemäß den gesetzlichen Bestimmungen den zuständigen Sammelpunkten zu. Endnutzer sind in Deutschland gesetzlich zur Rückgabe von Altbatterien an einer geeigneten Annahmestelle verpflichtet. Batterien können im Handelsgeschäft unentgeltlich zurückgegeben werden.



Indem Sie diese Produkte und Batterien ordnungsgemäß entsorgen, helfen Sie dabei, wertvolle Ressourcen zu schützen und eventuelle negative Auswirkungen auf die menschliche Gesundheit und die Umwelt zu vermeiden.

Für mehr Informationen zu Sammlung und Recycling, wenden Sie sich bitte an Ihren örtlichen Abfallentsorgungsdienstleister.

Gemäß Landesvorschriften können wegen nicht ordnungsgemäßer Entsorgung dieses Abfalls Strafgebühren verhängt werden.



### Hinweis für das Batteriesymbol (Symbol unten):

Dieses Symbol kann in Kombination mit einem chemischen Symbol abgebildet sein. In diesem Fall erfolgt dieses auf Grund der Anforderungen derjenigen Richtlinien, die für die betreffende Chemikalie erlassen wurden.

**Pb**



# DISPOSAL OF UNIT

(French)

## L'élimination des équipements et des batteries usagés

Applicable uniquement dans les pays membres de l'Union européenne et les pays disposant de systèmes de recyclage.



Apposé sur le produit lui-même, sur son emballage, ou figurant dans la documentation qui l'accompagne, ce pictogramme indique que les piles, appareils électriques et électroniques usagés, doivent être séparées des ordures ménagères.



Afin de permettre le traitement, la valorisation et le recyclage adéquats des piles et des appareils usagés, veuillez les porter à l'un des points de collecte prévus, conformément à la législation nationale en vigueur.



En les éliminant conformément à la réglementation en vigueur, vous contribuez à éviter le gaspillage de ressources précieuses ainsi qu'à protéger la santé humaine et l'environnement.

Pour de plus amples renseignements sur la collecte et le recyclage, veuillez vous renseigner auprès des collectivités locales.



Le non-respect de la réglementation relative à l'élimination des déchets est passible d'une peine d'amende.

### Note relative au pictogramme à apposer sur les piles (pictogramme du bas) :

Si ce pictogramme est combiné avec un symbole chimique, il répond également aux exigences posées par la Directive relative au produit chimique concerné.

**Pb**

(Spanish)

## Eliminación de Aparatos Viejos y de Pilas y Baterías

Solamente para la Unión Europea y países con sistemas de reciclado.



Estos símbolos en los productos, su embalaje o en los documentos que los acompañen significan que los productos eléctricos y electrónicos y pilas y baterías usadas no deben mezclarse con los residuos domésticos.



Para el adecuado tratamiento, recuperación y reciclaje de los productos viejos y pilas y baterías usadas llévelos a los puntos de recogida de acuerdo con su legislación nacional. En España, los usuarios están obligados a entregar las pilas en los correspondientes puntos de recogida. En cualquier caso, la entrega por los usuarios será sin coste alguno para éstos. El coste de la gestión medioambiental de los residuos de pilas, acumuladores y baterías está incluido en el precio de venta.



Si los elimina correctamente ayudará a preservar valuosos recursos y evitará potenciales efectos negativos sobre la salud de las personas y sobre el medio ambiente.

Para más información sobre la recogida u reciclaje, por favor contacte con su ayuntamiento.

Puede haber sanciones por una incorrecta eliminación de este residuo, de acuerdo con la legislación nacional.

### Nota para el símbolo de pilas y baterías (símbolo debajo):

Este símbolo puede usarse en combinación con el símbolo químico. En este caso, cumple con los requisitos de la Directiva del producto químico indicado.

**Pb**

# DISPOSAL OF UNIT

(Portuguese)

## Eliminação de Equipamentos Usados e Baterias

Apenas para a União Europeia e países com sistemas de reciclagem



Estes símbolos nos produtos, embalagens, e/ou documentos que os acompanham indicam que os produtos elétricos e eletrónicos e as baterias usados não podem ser misturados com os resíduos urbanos indiferenciados.



Para um tratamento adequado, reutilização e reciclagem de produtos e baterias usados, solicitamos que os coloque em pontos de recolha próprios, em conformidade com a respetiva legislação nacional.



Ao eliminar estes produtos corretamente estará a ajudar a poupar recursos valiosos e a prevenir quaisquer potenciais efeitos negativos sobre o ambiente e a saúde humana.

Para mais informações acerca da recolha e reciclagem, por favor contacte a sua autarquia local.



De acordo com a legislação nacional podem ser aplicadas contraordenações pela eliminação incorreta destes resíduos.

### Nota para o símbolo da bateria (símbolo na parte inferior):

Este símbolo pode ser utilizado conjuntamente com um símbolo químico. Neste caso estará em conformidade com o estabelecido na Diretiva referente aos produtos químicos em causa.

**Pb**

(Italian)

## Smaltimento di vecchie apparecchiature e batterie usate

Solo per Unione Europea e Nazioni con sistemi di raccolta e smaltimento



Questi simboli sui prodotti, sull'imballaggio e/o sulle documentazioni o manuali accompagnanti i prodotti indicano che i prodotti elettrici, elettronici e le batterie usate non devono essere buttati nei rifiuti domestici generici.



Per un trattamento adeguato, recupero e riciclaggio di vecchi prodotti e batterie usate vi invitiamo a portarli negli appositi punti di raccolta secondo la legislazione vigente nel vostro paese.



Con uno smaltimento corretto, contribuirete a salvare importanti risorse e ad evitare i potenziali effetti negativi sulla salute umana e sull'ambiente.

Per ulteriori informazioni su raccolta e riciclaggio, vi invitiamo a contattare il vostro comune.



Lo smaltimento non corretto di questi rifiuti potrebbe comportare sanzioni in accordo con la legislazione nazionale.

### Note per il simbolo batterie (simbolo sotto):

Questo simbolo può essere usato in combinazione con un simbolo chimico. In questo caso è conforme ai requisiti indicati dalla Direttiva per il prodotto chimico in questione.

**Pb**

# DISPOSAL OF UNIT

(Dutch)

## Het ontdoen van oude apparatuur en batterijen.

Enkel voor de Europese Unie en landen met recycle systemen.



Deze symbolen op de producten, verpakkingen en/of begeleidende documenten betekenen dat gebruikte elektrische en elektronische producten en batterijen niet samen mogen worden weggegooid met de rest van het huishoudelijk afval.



Voor een juiste verwerking, hergebruik en recycling van oude producten en batterijen, gelieve deze in te leveren bij de desbetreffende inleverpunten in overeenstemming met uw nationale wetgeving.



Door ze op de juiste wijze weg te gooien, helpt u mee met het besparen van kostbare hulpbronnen en voorkomt u potentiële negatieve effecten op de volksgezondheid en het milieu.

Voor meer informatie over inzameling en recycling kunt u contact opnemen met uw plaatselijke gemeente.



Afhankelijk van uw nationale wetgeving kunnen er boetes worden opgelegd bij het onjuist weggooien van dit soort afval.

**Pb**

**Let op: het batterij symbool (Onderstaand symbool).**

Dit symbool kan in combinatie met een chemisch symbool gebruikt worden. In dit geval volstaan de eisen, die zijn vastgesteld in de richtlijnen van de desbetreffende chemische stof.

(Swedish)

## Avfallshantering av produkter och batterier

Endast för Europeiska Unionen och länder med återvinningssystem



Dessa symboler på produkter, förpackningar och/eller medföljande dokument betyder att förbrukade elektriska och elektroniska produkter och batterier inte får blandas med vanliga hushållssopor.



För att gamla produkter och använda batterier ska hanteras och återvinnas på rätt sätt ska dom lämnas till passande uppsamlingsställe i enlighet med nationella bestämmelser.



Genom att ta göra det korrekt hjälper du till att spara värdefulla resurser och förhindrar eventuella negativa effekter på människors hälsa och på miljön.

För mer information om insamling och återvinning kontakta din kommun.



Olämplig avfallshantering kan beläggas med böter i enlighet med nationella bestämmelser.

**Pb**

**Notering till batterisymbolen (nedanför):**

Denna symbol kan användas i kombination med en kemisk symbol. I detta fall uppfyller den de krav som ställs i direktivet för den aktuella kemikalien.

# DISPOSAL OF BATTERY

## Location of a nickel-metal-hydride battery

This unit is provided a nickel-metal hydride battery for the power failure alarm. The battery is located in the control box inside the LED cover. (Fig. 1)




The high voltage components are enclosed in the control box. The cover should be removed by a qualified engineer or a service personnel only to prevent the electric shock.

## Disposal of nickel-metal-hydride battery

1. Turn off the power supply switch and disconnect the power supply plug.
2. As shown in the Fig. 1, remove 2 fixing screw A at the front lower of the LED cover with a screw driver.
3. Remove 2 screw B at the back of the foamed side table assy. (Fig. 2)
4. After disconnecting the connector, remove the LED cover. Then by unscrewing 2 screw C on the top of the foamed side table assy, remove the foamed side table assy. To remove the control box cover, unscrew 4 fixing screw D. (Fig. 3)
5. Remove 2 screw E fixing the battery mounting plate securing the battery and then disconnect the connector. (Fig. 4)
6. Take out the battery.

## Handling of battery

Cover the battery terminal with an insulating tape to avoid the short circuit. Then follow the procedure for recycling or proper disposal.




# TEMPERATURE RECORDER

## WARNING




Disconnect the power supply plug before attaching the temperature recorder or it may cause electric shock or fire.

The figure below shows the description of a temperature recorder.



### Feeding of chart


1. Open the door and let down the pen holder lever; the pen point is apart from the chart. (Fig. 1)
2. Pull the cartridge out of the mounted position. (Fig. 2)
3. Set a new chart in place on the rear bottom of the cartridge. Set the hole on the chart in the cog of the chart driving assembly and feed the chart in the direction of the arrow by driving the cog wheel.
4. Set the chart according to the mark of day and time. (Fig. 3)
5. When mounting, lay down the cartridge first, and push into the mounting position with the groove set a the projection on the position.
6. Set up the cartridge and settle in position.
7. Lift up the pen holder lever and close the door.



### Replacement of cell

Replace the dry cell once a year as follows:

1. Put down the pen holder lever first, then pull out cartridge from the mounted position.
2. Reset the battery box cover with wire, on the bottom left side, and take out the cell.
3. Set a new cell in the case.
4. Shut the battery box cover after exchanging cells. Mount the cartridge as before and lay up the penholder lever.



### Setting of backup temperature


1. There is a warning indicator (red guide) on the top of the recording indicator. Adjust the warning indicator with your finger to backup temperature, at which temperature the backup cooling kit starts to operate.
2. Set at a temperature 15°C higher than that in the freezer chamber.

# TEMPERATURE RECORDER

## Pen holder

Install a ink pen in the recording hand as shown in the figure. Make sure that the ink pen is completely inserted for accurate recording. The ink pen is packed together with the chart.

1. Let down the penholder lever, then pull out the cartridge from the mounted position.
2. Insert a new pen.
3. Mount the cartridge into place.
4. Lift up the penholder lever and make sure that the pen point touches the chart.



# BACKUP COOLING KIT

## Backup cooling kit for MDF-1156ATN (Liquid N<sub>2</sub> type)

1. Connect the liquid N<sub>2</sub> tank with the backup cooling kit joint by using connect tube and the joint attached with this freezer. This work should be done by high pressure gas works specialists. The pressure of liquid N<sub>2</sub> tank should be regulated at 0.5 kg/cm<sup>2</sup>G.
2. Operate the freezer until the freezer temperature reaches the required level.
3. Set the warning indicator of the temperature recorder at 15°C higher than the freezer temperature.
4. Switch on the liquid N<sub>2</sub> backup. This completes the backup cooling kit setting.
5. Make sure that liquid N<sub>2</sub> spouts into the freezer chamber with the backup test switch.
6. This completes the setting of the tank for automatic spouting in case the freezer temperature should match the alarm set temperature.

# SPECIFICATIONS

| Product name           | Ultra-Low Temperature Freezer<br>MDF-1156                                                                                                           | Ultra-Low Temperature Freezer<br>MDF-1156ATN                                                     |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| External dimensions    | W1400 mm x D800 mm x H945 mm                                                                                                                        |                                                                                                  |
| Internal dimensions    | W500 mm x D450 mm x H572 mm                                                                                                                         |                                                                                                  |
| Effective capacity     | 128 L                                                                                                                                               |                                                                                                  |
| Exterior               | Painted steel                                                                                                                                       |                                                                                                  |
| Interior               | Aluminum plate                                                                                                                                      |                                                                                                  |
| Door                   | Painted steel                                                                                                                                       |                                                                                                  |
| Insulation             | Rigid polyurethane foamed-in place                                                                                                                  |                                                                                                  |
| Compressor             | Hermetic type, 1100 W x 2                                                                                                                           |                                                                                                  |
| Condenser              | Fin and tube type (high stage side), Shell and tube type (low stage side)                                                                           |                                                                                                  |
| Evaporator             | Tube on sheet type                                                                                                                                  |                                                                                                  |
| Refrigerant            | R-407D (high stage side), HFC mixed refrigerant (low stage side)                                                                                    |                                                                                                  |
| Temperature controller | Microprocessor control system                                                                                                                       |                                                                                                  |
| Temperature display    | Digital display                                                                                                                                     |                                                                                                  |
| Sensor                 | Platinum resistance (Pt 100 Ω)                                                                                                                      |                                                                                                  |
| Alarm                  | buzzer, lamp, remote alarm contact                                                                                                                  |                                                                                                  |
| Battery                | For power failure alarm; Nickel-metal-hydride battery, DC 6 V, 650 mAh, Auto-recharge<br>For backup cooling kit; Lead storage battery, DC 6 V, 4 Ah |                                                                                                  |
| Temperature recorder   | ----                                                                                                                                                | Dry cell driven<br>Recording chart (for 2 months)                                                |
| Backup cooling kit     | ----                                                                                                                                                | Liquid N <sub>2</sub>                                                                            |
| Accessories            | 1 set of key, 2 rubber caps, 1 inner lid, 1 scraper                                                                                                 |                                                                                                  |
|                        | ----                                                                                                                                                | 6 recording paper, 1 dry cell<br>2 ink pens (cartridge)<br>1 connect tube for backup cooling kit |
| Weight                 | 265 kg                                                                                                                                              | 272 kg                                                                                           |
| Option                 | Storage case (MDF-49SC)                                                                                                                             |                                                                                                  |

**Note:** Design or specifications will be subject to change without notice.

The battery for power failure alarm is an article for consumption. It is recommended that the battery will be replaced about every 3 years. Contact our sales representative or agent at the time of replacement of the battery for recycling.

Refer to the updated catalog when ordering an optional component.

# PERFORMANCE

|                         |                                                                            |                                                                   |                               |
|-------------------------|----------------------------------------------------------------------------|-------------------------------------------------------------------|-------------------------------|
| Product name            | Ultra-Low Temperature Freezer<br>MDF-1156/MDF-1156ATN                      |                                                                   |                               |
| Model number            | MDF-1156-PB<br>MDF-1156ATN-PB                                              | MDF-1156-PK<br>MDF-1156ATN-PA<br>MDF-1156ATN-PK<br>MDF-1156ATN-PR | MDF-1156-PE<br>MDF-1156ATN-PE |
| Cooling performance     | -152°C at the center of the chamber; (Ambient temperature; +30°C, no load) |                                                                   |                               |
| Control range           | -130°C to -152°C                                                           |                                                                   |                               |
| Power voltage           | AC 220 V                                                                   | AC 220 V                                                          | AC 230 V/240 V                |
| Rated frequency         | 50 Hz                                                                      | 60 Hz                                                             | 50 Hz                         |
| Rated power consumption | 1550 W                                                                     | 1700 W                                                            | 1550 W/1600 W                 |
| Alarm duration          | 9 hours                                                                    |                                                                   |                               |
| Noise level             | 52 dB/54 dB (50 Hz/60 Hz) [A] (background noise; 20 dB)                    |                                                                   |                               |
| Maximum pressure        | 2844 kPa                                                                   |                                                                   |                               |

**Note:** The unit with CE mark complies with EU directives.



**⚠ CAUTION**

**Please fill in this form before servicing.  
Hand over this form to the service engineer to keep for his and your safety.**

## Safety check sheet

1. Freezer contents :

- Risk of infection: Yes No  
Risk of toxicity: Yes No  
Risk from radioactive sources: Yes No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit

Unit interior

- No contamination Yes No  
Decontaminated Yes No  
Contaminated Yes No

Others:

3. Instructions for safe repair/maintenance/disposal of the unit

- a) The unit is safe to work on Yes No  
b) There is some danger (see below) Yes No

Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

|                                                   |                |                |                       |
|---------------------------------------------------|----------------|----------------|-----------------------|
| Product name:<br>Ultra-low temperature<br>freezer | Model:<br>MDF- | Serial number: | Date of installation: |
|---------------------------------------------------|----------------|----------------|-----------------------|

Please decontaminate the unit yourself before calling the service engineer.





# PHC Corporation

1-1-1 Sakada, Oizumi-machi, Ora-gun, Gunma 370-0596, Japan

© PHC Corporation 2018

Printed in Japan  
LDCL027700-0  
S0418-0