## SUB Aqua Pro unstirred water bath range Accessories



Accessories										
	SAP2	SAP2S	SAP5	SAP12	SAP18	SAP26	SAP34	SAPD		
Capacity	2L	2L	5L	12L	18L	26L	34L	5L and 12L		
	Replacement polycarbonate transparent lids*									
	AQL2	AQL5	AQL5	AQL12	AQL26	AQL26	-	AQL5, AQL12		
	Directs condensation away from immersed vessels, avoids contamination, reduces evaporation and saves energy									
	Stainless stee	el sloping lids*								
	-	LU6	LU6	LU14	LU28	LU28	LU36	LU6, LU14		
(logar) (logar)	Flat lids*		. ==		. 500	1.500	1.500			
	-	-	LF6 (2 ring sets)	LF14 (4 ring sets)	LF28 (6 ring sets)	LF28 (6 ring sets)	LF36 (8 ring sets)	LF6 / LF14		
	With ring sets of variable hole diameter to accommodate tall vessels whilst reducing evaporation									
	Polypropylene	e spheres* (pa								
	1 x PS20	1 x PS20	1 x PS20	1 x PS20	2 x PS20	2 x PS20	3 x PS20	2 x PS20		
	Useful alternative to a lid, minimises evaporation and heat loss whilst allowing easy access to vessels in the bath; particularly useful for tall vessels									
	Raised shelves – reversible, allows two shelf depths. h = shelf height above tank base (mm)									
				RS14H (h 40 or 78)	RS18H (h 40 or 135)	RS28H (h 45 or 135)	RS36H h 45 or 135)	RS14H (h 40 or 78)		
	-	-	-	shelf covers	shelf covers	shelf covers	shelf covers	shelf covers		
				half area of SAP12	half area of SAP18	half area of SAP26	half area of SAP34	half area of SAPD		
	Racks** (no. p	per bath)								
	-	-	1 x J2	2 x J2	4 x J2	4 x J2	6 x J2	1 + 2 x J2		
	Choice of 8 variants to accommodate different tube diameters and microtubes (see below)									
	Replacement	base trays								
	AQBT2	AQBT5	AQBT5	AQBT12	AQBT26	AQBT26	SBT36	AQBT5 & AQBT12		
	Required if flat-bottomed flasks are to be placed directly on the base of the bath and to promote thermal convection in the bath									
-8	Drain stopper	•								
	-	-	-	UWB-DS (pack of 5)	UWB-DS (pack of 5)	UWB-DS (pack of 5)	UWB-DS (pack of 5)	UWB-DS (pack of 5)		
	Recommended when using heat transfer beads, to prevent beads entering drain hole									

 $<sup>^{\</sup>star}$  lid or spheres recommended for use above 60°C

Unstirred water bath racks**										
J2 racks	Tube size Ø	Capacity	J2 racks	Tube size Ø	Capacity					
J2-10	10 mm	84	J2-25	25 mm	18					
J2-13	13 mm	55	J2-30	30 mm	12					
J2-16	16 mm	36	J2-SE	0.5 ml	105					
J2-19	19 mm	32	J2-LE	1.5 ml	65					